

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com

www.ucytech.com

Data sheet Marmoset IFN-γ ELISA antibody pair; 10-plate format

Cat. No.: CT770-10

Coating antibodies (2 vials)

Product: Monoclonal antibody to marmoset interferon gamma (IFN-γ)

Isotype: Mouse IgG₁

Production: In vitro using serum free medium

Purification: DEAE ion exchange chromatography

Buffer: Prior to lyophilization: 0.25 ml PBS + 125 mM trehalose

Application: Inject 0.25 ml distilled water into the vial and dilute 100 times in PBS. The content

of one vial is sufficient for five 96-well ELISA plates (480 determinations; 50

μl/well).

Detection antibodies (2 vials)

Product: Biotinylated monoclonal antibody to marmoset interferon gamma (IFN-γ)

Isotype: Mouse IgG₁

Production: In vitro using serum free medium

Purification: DEAE ion exchange chromatography

Labeling: With Biotin-7-NHS (N-hydroxysuccinimide)

Buffer: Prior to lyophilization: 0.5 ml PBS + 1% BSA + 125 mM trehalose

Application: Inject 0.5 ml distilled water into the vial and dilute 100 times in PBS + 0.5% BSA +

0.05% Tween-20. The content of one vial is sufficient for five 96-well ELISA plates

(480 determinations; 100 μl/well).

Standards (5 vials)

Product: Recombinant marmoset interferon gamma (IFN-γ)

Application: Cytokine standard in an ELISA system

Reconstitution: Inject 0.5 ml distilled water into the vial. Use immediately.

Conjugate (2 vials)

Product: SPP conjugate (Streptavidin-HRP)

Application: Inject 0.5 ml distilled water into the vial and dilute 100 times in PBS + 0.5% BSA +

0.05% Tween-20. The content of one vial is sufficient for five 96-well ELISA plates

(480 determinations; 100 µl/well).

The product should be used in combination with TMB substrate.

General

Sensitivity: 1 pg/ml

Specificity: Validated for detecting natural and recombinant marmoset IFN-γ. Standards are

calibrated with recombinant marmoset IFN-y.

Sterility: Membrane filtered (0.2 µm)

Stability: Lyophilized SPP conjugate is stable for at least one year at -20°C in the dark, the

other lyophilized products are stable for more than one year at $4\,^{\circ}\text{C}$.

After reconstitution, the antibodies are stable for one year at $4\,^\circ\text{C}$ (if kept sterile) and SPP for minimal one year at $\text{-}20\,^\circ\text{C}$ in the dark. The reconstituted standard

preparation should be used immediately.



For research use only Version 140820