


U-CyTech B.V.
Yalelaan 48
3584 CM Utrecht
The Netherlands
+31 85 073 1460
info@ucytech.com
www.ucytech.com

Data sheet Human IFN- γ /IL-2 FluoroSpot kit

Catalogue no. CT852-F2 (2-plate format)

Species reactivity human



Typical example of IFN- γ /IL-2 specific spots produced by 2.5×10^3 human PBMC
In green: IFN- γ , in red: IL-2, in yellow: IFN- γ /IL-2.

Preincubation: 42h
Incubation FluoroSpot: 20h
Stimulation: PMA/ionomycin

References can be found in the Reference Database at www.ucytech.com/references.

