


U-CyTech B.V.

Yalelaan 42 3584 CM Utrecht The Netherlands +31 85 073 1460 info@ucytech.com www.ucytech.com Data sheet Human IFN-γ/IL-5 FluoroSpot kit

Catalogue no. CT853-F2 (2-plate format)

Species reactivity human



Typical example of IFN- γ /IL-5 specific spots produced by 2.5x10⁴ human PBMC In green: IFN- γ , in red: IL-5, in yellow: IFN- γ /IL-5.

Preincubation: 42 h Incubation FluoroSpot: 20 h

Stimulation: PMA/ionomycin