


## U-CyTech B.V.

Yalelaan 42 3584 CM Utrecht The Netherlands +31 85 073 1460 info@ucytech.com www.ucytech.com Data sheet Human IFN-γ/IL-1β FluoroSpot kit

Catalogue no. CT856-F2 (2-plate format)

Species reactivity: human



Typical example of IFN-y/IL-1B specific spots produced by  $5x10^4$  human PBMC: In green: IFN-y, in red: IL-1B, in yellow: IFN-y/IL-1B.

Preincubation: 42 h Incubation FluoroSpot: 20 h

Stimulation: PMA/ionomycin