


U-CyTech B.V.
Yalelaan 42
3584 CM Utrecht
The Netherlands
+31 85 073 1460
info@ucytech.com
www.ucytech.com

Data sheet Human IFN- γ /IL-4 FluoroSpot kit
Catalogue no. CT857-F2 (2-plate format)

Species reactivity human



Typical example of IFN- γ /IL-4 specific spots produced by 5×10^4 human PBMC
In green: IFN- γ , in red: IL-4, in yellow: IFN- γ /IL-4.

Preincubation: 42 h
Incubation FluoroSpot: 20 h
Stimulation: PMA/ionomycin

