

Data sheet Mouse IFN- γ ELISPOT antibody pair; 20-plate format

Cat. No.: CT655-20

Coating antibodies (1 vial)

Product: Monoclonal antibody to mouse interferon gamma (IFN- γ)
Isotype: Mouse IgG₁
Production: *In vitro* using serum free medium
Purification: Ammonium sulphate precipitation and DEAE ion exchange chromatography
Contents: Each vial contains sufficient material for coating of twenty 96-well ELISPOT plates
Buffer: Prior to lyophilization: 1.0 ml PBS + 125 mM trehalose
Application: Coating antibody in an ELISPOT system
Reconstitution: Dissolve the contents of the vial by injection of 1.0 ml distilled water into the vial and dilute 100 times in PBS. The total amount of one vial is sufficient for twenty 96-well ELISPOT plates (1920 determinations; 50 μ l/well).

Detection antibodies (1 vial)

Product: Biotinylated polyclonal antibody to mouse interferon gamma (IFN- γ)
Isotype: Rabbit IgG
Purification: Ammonium sulphate precipitation, protein A- and ligand-affinity chromatography
Labeling: With Biotin-7-NHS (N-hydroxysuccinimide)
Contents: Each vial contains sufficient material for twenty 96-well ELISPOT plates
Buffer: Prior to lyophilization: 2.0 ml PBS + 1% BSA + 125 mM trehalose
Application: Detection antibody in an ELISPOT system
Reconstitution: Dissolve the contents of the vial by injection of 2.0 ml distilled water into the vial and dilute 100 times in Dilution buffer (see Technical Data Sheet). The total amount of one vial is sufficient for twenty 96-well ELISPOT plates (1920 determinations; 100 μ l/well).

General

Specificity: Validated for detecting mouse IFN- γ
Sterility: Membrane filtered (0.2 μ m)
Stability: The lyophilized products are stable for at least one year at 4°C (expiry date is indicated on the vials).
After reconstitution, the antibodies are stable for several months at 4°C (if kept sterile) or for minimal one year at -20°C.
References: Engler, O.B. *et al.* 2004 Vaccine 23: 58-68
Hanke, T. *et al.* 2003. J. Gen. Virol. 84: 361-368
Hu, H. *et al.* 2007. Clin. Vaccine Immunol. 14: 894-901
Jin, Y. *et al.* 2005. Clin. Diagn. Lab. Immunol. 12: 1364-1369
Li, S. *et al.* 2011. Vaccine 29: 3695-3702