

Data sheet Mouse IL-5 ELISPOT antibody pair; 10-plate format

U-CyTech BV Yalelaan 48 3584 CM Utrecht The Netherlands P +31.30.253 5960 F +31.30.253 9344 INFO@ucytech.com

| Cat. No.: | CT660-10 | |
|------------------------------|--|--|
| Coating antibodies (2 vials) | | |
| Product: | Monoclonal antibody to mouse interleukin 5 (IL-5) | |
| Isotype: | Rat IgG ₁ | |
| Production: | In vitro using serum free medium | |
| Purification: | Ion exchange chromatography | |
| Contents: | Each vial contains sufficient material for coating of five 96-well ELISPOT plates | |
| Buffer: | Prior to lyophilization: 0.25 ml PBS + 125 mM trehalose | |
| Application: | Coating antibody in an ELISPOT system | |
| Reconstitution: | Dissolve the contents of one vial by injection of 0.25 ml distilled water into the vial and dilute 100 times in PBS. The total amount of one vial is sufficient for five 96-well ELISPOT plates (480 determinations; 50 μ l/well). | |

Detection antibodies (2 vials)

| Product: | Biotinylated monoclonal antibody to mouse interleukin 5 (IL-5) |
|-----------------|--|
| Isotype: | Rat IgG ₁ |
| Production: | In vitro using serum free medium |
| Purification: | Protein G-affinity chromatography |
| Labeling: | With Biotin-7-NHS (N-hydroxysuccinimide) |
| Contents: | Each vial contains sufficient material for five 96-well ELISPOT plates |
| Buffer: | Prior to Iyophilization: 0.5 ml PBS + 1% BSA + 125 mM trehalose |
| Application: | Detection antibody in an ELISPOT system |
| Reconstitution: | Dissolve the contents of one vial by injection of 0.5 ml distilled water into the vial and dilute 100 times in Dilution buffer (see Technical Data Sheet). The total amount of one vial is sufficient for five 96-well ELISPOT plates (480 determinations; 100μ /well). |

General

| Specificity: | Validated for detecting mouse IL-5 |
|--------------|--|
| Sterility: | Membrane filtered (0.2 µm) |
| Stability: | The lyophilized products are stable for at least one year at 4°C (expiry date is indicated on the vials). After reconstitution, the antibodies are stable for several months at 4°C (if kept sterile) or for minimal one year at -20°C. |

