

Data sheet Human Perforin ELISPOT antibody pair; 10-plate format

U-CyTech BV Yalelaan 48 3584 CM Utrecht The Netherlands P +31.30.253 5960 F +31.30.253 9344 INFO@ucytech.com www.ucytech.com

Cat. No.:	CT653-10	
Coating antibodies (2 vials)		
Product:	Monoclonal antibody to human perforin	
lsotype:	Mouse IgG ₁	
Production:	In vitro using serum free medium	
Purification:	DEAE ion exchange chromatography	
Contents:	Each vial contains sufficient material for coating of five 96-well ELISPOT plates	
Buffer:	Prior to lyophilization: 0.25 ml PBS + 125 mM trehalose	
Application:	Coating antibody in an ELISPOT system	
Reconstitution:	Dissolve the contents of the vial by injection of 0.25 ml distilled water into the vial and dilute 100 times in PBS. The total amount of one vial is sufficient for five 96-well ELISPOT plates (480 determinations; 50 μ l/well).	

Detection antibodies (2 vials)

Product:	Biotinylated monoclonal antibody to human perforin
Isotype:	Mouse IgG ₁
Production:	In vitro using serum free medium
Purification:	DEAE ion exchange chromatography
Labeling:	With Biotin-7-NHS (N-hydroxysuccinimide)
Contents:	Each vial contains sufficient material for five 96-well ELISPOT plates
Buffer:	Prior to lyophilization: 0.5 ml PBS + 1% BSA + 125 mM trehalose
Application:	Detection antibody in an ELISPOT system
Reconstitution:	Dissolve the contents of the vial by injection of 0.5 ml distilled water into the vial and dilute 100 times in Dilution buffer (see Technical Data Sheet). The total amount of one vial is sufficient for twenty 96-well ELISPOT plates (480 determinations; 100μ /well).

General

Specificity:	Validated for detecting human perforin
Sterility:	Membrane filtered (0.2 µm)
Stability:	The lyophilized products are stable for at least one year at 4°C (expiry date is indicated on the vials). After reconstitution, the antibodies are stable for several months at 4°C (if kept sterile) or for minimal one year at -20°C.

