


U-CyTech B.V.
Yalelaan 48
3584 CM Utrecht
The Netherlands
+31 85 073 1460
info@ucytech.com
www.ucytech.com


Data sheet Human IgG B cell ELISPOT kit

Catalogue no. CT780-PR5 (5-plate format)

Species reactivity human



Total number of IgG secreting cells
 2×10^3 cells/well (Assay I)



Tetanus toxoid specific IgG secreting cells
 1×10^5 cells/well (Assay II)

Example of IgG specific spots produced by human PBMC

Individual was intramuscularly vaccinated with tetanus toxoid more than 15 years ago.

B cell activation: 5 days preincubation (2×10^6 PBMC/ml) in the presence of recombinant IL-2 and R848.

Incubation ELISPOT: 20 hours

References can be found in the Reference Database at www.ucytech.com/references.

