


U-CyTech B.V.
Yalelaan 42
3584 CM Utrecht
The Netherlands
+31 85 073 1460
info@ucytech.com
www.ucytech.com

Data sheet Monkey IgM B cell ELISPOT kit

Catalogue no. CT786-B5 (5-plate format)

Species reactivity rhesus macaque



Total number of **spontaneously** IgM secreting cells (5×10^4 cells/well)

Total number of IgM secreting cells after R848/IL-2 (5×10^3 cells/well) 5×10^4 cells/well

2.5×10^4 cell/well

Example of IgM specific spots produced by rhesus macaque PBMC

B cell activation: 5 days preincubation (2×10^6 PBMC/ml) in the absence (left) or presence (right) of recombinant IL-2 and R848.
Incubation on ELISPOT plate: 24 hours

