


U-CyTech B.V.
Yalelaan 42
3584 CM Utrecht
The Netherlands
+31 85 073 1460
info@ucytech.com
www.ucytech.com

Data sheet	Monkey IFN- γ /TNF- α FluoroSpot kit
Catalogue no.	CT901-F2 (2-plate format)
Species reactivity	rhesus macaque, cynomolgus monkey, pig-tailed macaque, baboon, African green monkey



Typical example of IFN- γ /TNF- α specific spots produced by 2.5×10^4 rhesus monkey PBMC
In green: IFN- γ , in red: TNF- α , in yellow: IFN- γ /TNF- α .

Preincubation:	42 h
Incubation FluoroSpot:	20 h
Stimulation:	PMA/ionomycin

